

# MARISORB®

## Qualities

- Marisorb® is a synthetic, monofilament, absorbable suture made from Polydioxanon.
- The tensile strength of Marisorb® is reduced by 50 per cent after 6 to 7 weeks, after about 180 to 210 days degradation by hydrolysis is completed.

## **Application:**

- . Gastrointestinal
- . Peritoneum
- . Mucosa
- . Ligatures
- . Fascial
- . Intra-dermal
- . Muscles
- . Subcutaneous
- . Arterial and Cardiovascular
- . Nerves
- . Urological
- . Gynaecological
- . Tendons
- . Plastic










## How supplied

- needled sutures
- loops
- ligatures (upon request)

# MARISORB®

synthetic, absorbable  
monofilament









## Needled sutures

Needle description	Article-No.	Colour	USP	metric	Length	UoS
 HR 13	17301035	violet	6/0	0,7	70 cm	24
	17301040	violet	5/0	1	70 cm	24
	17301044	violet	4/0	1,5	70 cm	24
	17301048	violet	3/0	2	70 cm	24
	17301052	violet	2/0	3	70 cm	24
 2xHR 13	17301094	violet	5/0	1	70 cm	24
 HR 17	17301174	violet	6/0	0,7	70 cm	24
	17301179	violet	5/0	1	70 cm	24
	17301183	violet	4/0	1,5	70 cm	24
	17301187	violet	3/0	2	70 cm	24
 HR 17s	17301199	violet	3/0	2	70 cm	24
	17301201	violet	2/0	3	70 cm	24
 2xHR 17	17301229	violet	3/0	2	90 cm	24
	17301232	violet	2/0	3	90 cm	24
	17301219	violet	4/0	1,5	70 cm	24
 HR 22	17301258	violet	4/0	1,5	70 cm	24
	17301299	violet	4/0	1,5	90 cm	24
	17301260	violet	3/0	2	70 cm	24
	17301264	violet	2/0	3	70 cm	24
	17301270	violet	0	3,5	70 cm	24
 HR 22s	17301355	violet	2/0	3	70 cm	24
	17301357	violet	0	3,5	70 cm	24
 HR 26	17301370	violet	4/0	1,5	70 cm	24
	17301372	violet	3/0	2	70 cm	24
	17301374	violet	2/0	3	70 cm	24
	17301376	violet	0	3,5	70 cm	24
 2xHR 26	17301400	violet	4/0	1,5	70 cm	24
	17301402	violet	3/0	2	70 cm	24
	17301404	violet	2/0	3	70 cm	24
	17301406	violet	0	3,5	70 cm	24

# MARISORB®

synthetic, absorbable  
monofilament












## Needled sutures

Needle description	Article-No.	Colour	USP	metric	Length	UoS
 HR 27s	17301540	violet	2/0	3	70 cm	24
	17301542	violet	0	3,5	70 cm	24
	17301545	violet	1	4	70 cm	24
 HR 30	17301600	violet	3/0	2	70 cm	24
 HR 37s	17301656	violet	2/0	3	70 cm	24
	17301670	violet	0	3,5	70 cm	24
	17301672	violet	1	4	70 cm	24
	17301701	violet	1	4	90 cm	24
 HR 43	17301772	violet	2/0	3	90 cm	24
	17301774	violet	0	3,5	90 cm	24
	17301760	violet	1	4	70 cm	24
	17301776	violet	1	4	90 cm	24
 HR 43s	17301790	violet	0	3,5	90 cm	24
	17301793	violet	1	4	90 cm	24
 HR 48s	17301931	violet	1	4	90 cm	24
	17301933	violet	2	5	90 cm	24
 HRT 17	17302200	violet	5/0	1	70 cm	24
	17302204	violet	4/0	1,5	70 cm	24
	17302205	violet	3/0	2	70 cm	24
	17302209	violet	2/0	3	70 cm	24
 HRT 26	17302310	violet	2/0	3	70 cm	24

# MARISORB®

synthetic, absorbable  
monofilament











## Needed sutures

Needle description	Article-No.	Colour	USP	metric	Length	UoS
 2xDR 10	17302540	violet	6/0	0,7	70 cm	24
 DR 12	17302632 17302644 17302645	violet violet violet	6/0 6/0 5/0	0,7 0,7 1	45 cm 70 cm 70 cm	24 24 24
 2xDR 12	17302658 17302660	violet violet	6/0 4/0	0,7 1,5	70 cm 70 cm	24 24
 DS 12	17302881 17302882	violet violet	5/0 4/0	1 1,5	45 cm 45 cm	24 24
 2xDS 12	17302914 17302916	violet violet	6/0 5/0	0,7 1	70 cm 70 cm	24 24
 DS 16	17303002 17303004	violet violet	5/0 3/0	1 2	45 cm 45 cm	24 24
 DS 19	17303035 17303037 17303055 17303059	violet violet violet violet	4/0 3/0 4/0 3/0	1,5 2 1,5 2	45 cm 45 cm 70 cm 70 cm	24 24 24 24
 DS 24	17303180 17303184	violet violet	4/0 3/0	1,5 2	70 cm 70 cm	24 24
 DSM 13	17303312 17303314 17303316	violet violet violet	5/0 4/0 3/0	1 1,5 2	45 cm 45 cm 45 cm	24 24 24
 DSM 16	17303421 17303422	violet violet	5/0 4/0	1 1,5	45 cm 45 cm	24 24
 DSM 19	17303529 17303535	violet violet	4/0 4/0	1,5 1,5	45 cm 70 cm	24 24

# MARISORB®

synthetic, absorbable  
monofilament





Needled sutures

Needle description	Article-No.	Colour	USP	metric	Length	UoS
 DSM 24	17303570	violet	2/0	3	70 cm	24
 HS 12	17303668	violet	5/0	1	45 cm	24
 HS 30	17303711 17303713	violet violet	0 1	3,5 4	70 cm 70 cm	24 24
 HS 37	17303751	violet	1	4	70 cm	24
 HS 40	17303779	violet	1	4	70 cm	24
 2xGR 19	17303825	violet	4/0	1,5	45 cm	24
 2xGR 65	17303849 17303852	violet violet	0 1	3,5 4	70 cm 70 cm	24 24
 FR 27	17303880	violet	0	3,5	70 cm	24
 2xFR 27	17303910 17303914	violet violet	3/0 2/0	2 3	70 cm 90 cm	24 24
 2xFR 36	17303950	violet	2/0	3	70 cm	24

# MARISORB®

synthetic, absorbable  
monofilament

## Loop sutures

Needle description	Article-No.	Colour	USP	metric	Length	UoS
 HR 43s	17304222	violet	0	3,5	150 cm	12
	17304224	violet	1	4	150 cm	12
 HR 48s	17304335	violet	1	4	90 cm	12
	17304340	violet	1	4	150 cm	12
	17304351	violet	1	4	180 cm	12
 HRT 43	17304441	violet	1	4	150 cm	12
 HRT 48s	17304462	violet	1	4	150 cm	12

## Pretied ligating loop

	Article-No.	Colour	USP	metric		UoS
	17304781	violet	0	3,5		12
	17304782	violet	1	4		12